



ClearSkin PRO

Focus. Flexibility. Freedom. **That's Harmony.**





ClearSkin PRO

When it comes to rejuvenation, always work with a **PRO**.

ClearSkin PRO is one of Alma's breakthrough applicators, offering outstanding proven results for skin rejuvenation treatments through extraordinarily powerful energy emission. It is an innovative applicator, combining a non-ablative mid infra-red laser with simultaneous contact cooling to promote the synthesis of new skin tissue while maintaining a high level of safety. ClearSkin PRO is a significant evolutionary development that enables enhanced clinical capabilities alongside existing solutions, adding real value to your clinical offer.

With such a powerful technology in hand, it's no wonder the results are so amazing.

Powerful Technology for Positive Results

With **twice the power** of previous solutions, ClearSkin PRO applicator delivers up to 3000 mJ (61.2 mJ/Pixel) of non-ablative laser energy per pulse. This considerable increase in power provides a very significant clinical benefit and a proven, visible difference in skin rejuvenation treatments.


ClearSkin PRO's non-ablative Er:Glass 1540nm laser deeply penetrates the skin, creating columns of coagulation, which promotes new dermis synthesis and regenerates the extra cellular matrix while leaving the epidermis intact. The result is visibly younger looking skin.

ClearSkin PRO is launched on the new Harmony XL PRO Special Edition – a modern and refreshed exterior representing Alma's new perception of skin rejuvenation. It combines the strongest thermal effect with an already clinically proven array of applicators for holistic rejuvenation.

I have quite a lot of experience with various technologies and applicators, but ClearSkin PRO definitely offers something unique. Its high level of energy enables much more powerful treatments, and the results are clearly apparent. This is a significant development which adds real value to my clinic, and exceedingly better results for my patients.



A coagulation column of thermal damage, around 850 microns deep, is shown in the dermis while the epidermis remains intact (dark purple outer border). The surrounding unaffected dermal tissue will contribute to the synthesis of new collagen.







Courtesy of Alma's clinical department





Powerful Technology. Positive Results. Now, that's a PRO.

High power (up to 3000 mJ)

Proven and remarkable results

Zero downtime

Year-round treatments

Optimal & safe for all skin types

Exclusive Distributor:



PT. RADIAN ELOK DISTRIVERSA

Graha SKG-GN, Jl. Raya Kelapa Nias Blok GN-5, Kelapa Gading Barat, Jakarta Utara 14240, INDONESIA

T: +6221 45850088 www.radian-elok-distriversa.com

Alma Lasers GmbH

Nordostpark 100-102 90411 Nuremberg, Germany Tel. + 49 911 / 89 11 29-0 Fax + 49 911 / 89 11 29-99

Email: info@almalasers.com

www.almalasers.com

© 2020 Alma. All rights reserved. Alma ,its logo, Harmony, ClearSkin are trademarks and/or registered trademarks of Alma in the United States and/or other countries.

PBHA21111911 01

CONNECT WITH RED











